

IKA

designed for scientists

imLab
SCIENTIFIC EQUIPMENT
FOR LABORATORY & INDUSTRY

@ **Web contact**
www.imlab.eu
info@imlab.eu

☎ **Phone**
🇫🇷 +33(0)3 20 55 19 11
🇪🇸 +32(0)16 73 55 72

EN



IKA Specials Q1 / 2026
Innovative Laboratory Solutions

Take advantage of our time-limited offers on essential lab equipment.

Get your discounts now using the promo code **QC12026**

The discounts are valid until March 31st, 2026.



Patented, removable
shaking platform

Cooled Incubator and Incubator Shakers INC 125 FS & FC digital

- › Heating: up to 80 °C, cooling: down to -5 °C
- › Fan Speed: adjustable (0 – 100%) for optimal temperature distribution
- › Decontamination: 100 °C for 2.5 h
- › Cooling Incubator with eco-friendly refrigerant R290

2+1
Year
Warranty*



INC 125 FS digital (SP20)

Ident. No.: 0020117221

~~€ 5,666.00~~ **€ 4,816.10**

INC 125 FC digital Incubator cooled

Ident. No.: 0020126186

~~€ 4,944.00~~ **€ 4,202.40**

INC 125 FS digital (SP25)

Ident. No.: 0020117230

~~€ 5,666.00~~ **€ 4,816.10**

-15%

LIMITED Pipette Deal

**Every IKA PETTE vario 100 – 1000 µl
comes with a free tip box**

IKA PETTE vario 100 – 1000 µl

Ident. No.: 0020011216

€ 275.00



Free

NEW!


Handheld Dispenser T 10 battery ULTRA-TURRAX®

2+1
Year
Warranty*

- › Cordless, battery-operated homogenizer, ideal for mobile use
- › Volume Range: up to 100 ml (H₂O), fits microcentrifuge tubes
- › For on-site sampling, biosafety/laminar flow cabinets, tissue homogenisation
- › Brushless DC motor for quiet operation
- › Wide variety of dispersing tools available: blunt, sawtooth, stainless steel, plastic

Ident. No.: 0020103418

~~€ 1,112.00~~ **€ 889.60**
-20%

Best value in
its class


Rotary Evaporator RV 3 V Starter Solution

2+3
Year
Warranty*

Cost-effective complete setup includes:

- › HB eco heating bath (4 L, RT – 99 °C)
- › Vertical glass condenser (1 L flasks, 1500 cm² surface)
- › VACSTAR lite vacuum pump (min. 9 mbar)
- › VC 10 lite controller (manual & program modes)
- › RC 2 lite chiller (400 W, down to –10 °C, 0.35 bar pump)

ECO
energy-efficient
technology

Ident. No.: 0010013945

~~€ 7,632.00~~ **€ 6,105.60**
-20%

*coated V-C version available on request – please contact your IKA representative

With eco-friendly
refrigerant R 290


Compact chiller / circulator for cooling and heating RC 2 lite & HRC 2 lite

2+1
Year
Warranty*

- › **RC 2 lite:** 400 W cooling, ideal for evaporators, condensers, heat exchangers
- › **HRC 2 lite:** 400 W cooling + 1000 W heating, perfect for bioreactors, lab reactors, open baths
- › Fast temperature control with 1 L fill volume, 2.5 L working volume
- › Powerful pump (0.35 bar) supports multiple or open-bath applications

ECO
energy-efficient
technology

Combines best as peripheral devices:
EasySyn & HABITAT

RC 2 lite

Ident. No.: 0025006624

~~€ 2,736.00~~ **€ 2,323.60**
HRC 2 lite

Ident. No.: 0020104311

~~€ 3,136.00~~ **€ 2,665.60**
-15%

Bestseller



Magnetic Stirrer RCT basic ready-to-go Solution

- › Complete setup for precise temperature control
- › Includes stand rod, extension arm & boss head clamp for safe PT 1000 sensor positioning
- › Sensor placement optimized: outside stirring vortex, ≥ 20 mm immersion depth
- › Silicone cover (600 & 1000 ml beakers): splash guard, energy-saving lid, noise reduction

Ident. No.: 0010011525

~~€ 1,058.00~~ **€ 889.30**

-15%

2+1
Year
Warranty*

Compact &
user-friendly



Vacuum Pump VACSTAR lite Starter Solution

- › Includes VACSTAR lite pump + VC 10 lite controller
- › Integrated device communication with circulators for optimized workflows
- › Valve cleaning function extends system life & reduces maintenance

Ident. No.: 0010014496

~~€ 2,458.50~~ **€ 1,966.80**

-20%

2+1
Year
Warranty*

Ready-to-connect
packages available
for easy setup



Immersion and bridge circulators ICC 150 lite

- › High-performance temperature control unit for liquids up to 150 °C
- › Internal & external applications via integrated pressure/suction pump
- › Safety class III (DIN 12876), approved for flammable liquids
- › Overflow detection with automatic heating shutoff
- › Integrated stand for easy cleaning & component protection

Ident. No.: 0020105581

~~€ 969.00~~ **€ 775.20**

-20%

2+1
Year
Warranty*



Overhead Stirrer MINISTAR 40 digital

- › "Mini" Overhead stirrer for up to 25 L (H₂O)
- › Speed range: 50 – 1000 rpm with constant torque (max. 40 Ncm)
- › Handles viscosities up to 30,000 mPas
- › Protection class: IP54, suitable for ambient temperatures 5 – 40 °C

Ident. No.: 0025004886

~~€ 1,309.00~~ **€ 1,112.65**

-15%

2+8
Year
Warranty*

Space-saving design
with patented
temperature sensing



Disperser T 25 mini control Solution

- › High-performance disperser for up to 1.5 L (H₂O), 600 – 25,000 rpm
- › Ideal for low-viscosity sample homogenization in tight lab spaces
- › Digital speed control & brushless motor for quiet operation
- › Included accessories: dispersing tool, plate stand, boss head clamp

Ident. No.: 0010015125

~~€ 4,070.00~~ **€ 3,256.00**

-20%

2+1
Year
Warranty*

Free trial: 30-day full
access to all features &
700+ devices, no obligation



Laboratory Software labworldsoft® 6 Lite

- › Entry-level software for monitoring, controlling & documenting one IKA device
- › Ideal for automating lab processes & real-time calculations via PC
- › Supports FDA 21 CFR Part 11 compliance & advanced user management
- › Includes 1-year updates

Ident. No.: 0020017364

~~€ 1,184.00~~ **€ 888.00**

-25%