

IKA

designed for scientists

EN

Trusted Solutions



for Your Laboratory

IKA Specials Q2/ 2026
Innovative Laboratory Solutions

imLab



www.imlab.eu - info@imlab.eu



+33(0)3 20 55 19 11



+32(0)16 73 55 72

Take advantage of our time-limited offers on essential lab equipment.
The discounts are valid until June 30th, 2026.



reddot winner 2022



Universal and compact mixing and shaking device with magnetic drive

IKA MIXSTER SET 4

2+1
Year
Warranty*

- › Multi-position magnetic stirrer set with four drive units
- › For mixing up to 5 liters at a max speed of 3000 rpm per single unit
- › Can be expanded to up to 30 stirring positions
- › Can be transformed into a magnetic stirrer with 9 stirring positions for vessels with small volumes

Ident. No.: 0020027390

~~€ 1,620.00~~ **€ 1.458,00**

-10%

Shaker

MS 3 control

2+1
Year
Warranty*



- › Compact, universal small shaker suitable for shaking tasks with all small vessels and microtiter plates
- › Wide range of attachments
- › Attachment detection
- › Timer with countdown function
- › Stable in all speed ranges
- › Continuous or touch operation (with standard attachment)

Ident. No.: 0020016017

~~€ 1,091.00~~ **€ 981,90**

-10%

Mix, heat, cool –
all in one



+ Free
Pipette



Thermoshaker MATRIX Orbital Delta Plus

- › Large and clear display with a convenient menu structure
- › Firm quick-release fastener for quick and easy changing of the inserts/attachments
- › Quick and effective mixing of liquids with up to 3,000 rpm and 3 mm shaker diameter
- › Automatic performance of mixing, heating and cooling tasks through programmable ramps

Ident. No.: 0010006853

~~€ 2,839.00~~ € 2.555,10

-10%

2+1
Year
Warranty*



Mills MultiDrive basic

- › For a wide range of crushing tasks involving coarse and fine crushing
- › Excellent crushing performance with variable rotational speeds (from 3000 rpm to 20 000 rpm)
- › With interval operation at the press of a button
- › With integrated cooling for indirect heat dissipation
- › Variety of vessels available as accessories

Ident. No.: 0025002638

~~€ 3,835.00~~ € 3.451.50

-10%

2+1
Year
Warranty*

* Vessels are not part of the scope of delivery.



Single-Use Crushing ULTRA-TURRAX® Tube Drive P control Workstation

- › Powerful disperser for best stirring, dispersing, grinding and homogenisation
- › In particular for the safe processing of infectious, toxic or odour-intensive samples
- › With sealable and disposable tubes to eliminate the risk of cross-contamination
- › Programmable test conditions

Ident. No.: 00250085836

~~€ 1,820.00~~ € 1.638,00

-10%

2+1
Year
Warranty*

NEW!



Disperser Solution T 10 battery Complete Solution

2+1
Year
Warranty*

- › Complete solution for benchtop dispersing operation
- › Includes a reusable stainless steel dispersing tool (5 mm), a stand and a clamp
- › Large LCD display for easy parameter viewing
- › Brushless DC motor for quieter operation
- › Precise speed setting with numerical values
- › Timer and counter functions
- › Rechargeable battery with two charging options

Ident. No.: 0010015688

~~€ 2,335.00~~ **€ 1.868,00**

-20%



Disperser Solution T 18 digital Solution

2+1
Year
Warranty*

Complete Solution with T 18 digital disperser

Scope of delivery

- › T 18 digital ULTRA-TURRAX®
- › S 18 N - 19 G Dispersing tool
- › R 1826 Plate stand
- › R 182 Boss head clamp

Ident. No.: 0009024700

~~€ 2,549.00~~ **€ 2.294,10**

-10%

For stirring quantities
of up to 20 l (H₂O)



Overhead Stirrer RW 20 digital Package

2+1
Year
Warranty*

- › With digital display
- › Robust, slimline, ergonomic design
- › With constant power-drive
- › Two speed ranges for universal use from 60 – 2000 rpm
- › Push-through agitator shafts (only when stationary)

Ident. No.: 0009011600

~~€ 1,599.00~~ **€ 1.439,10**

-10%

Fully automated calibration by robots



Liquid Handling IKA PETTE Classic Kit

- › Mechanical air-cushion pipettes manufactured in Germany
- › Covers a volume range of 2 – 1000 µl
- › Fully autoclavable

Package consists of:

- › three single channel pipettes with variable volumes
- › Linea stand
- › two pipette tip boxes

Ident. No.: 0010013011

~~€ 713,00~~ € 641,70

-10%

2+1
Year
Warranty*

With natural refrigerant R290



Cooling and Heating Circulator HRC 2 GREEN basic

- › Speed-controlled refrigeration system with 620 W cooling capacity
- › With a working temperature range from -30°C to 100°C
- › With a powerful, continuously variable pressure and suction pump made of PEEK: (0.5 bar; 22 l/min)
- › 1500 W heating capacity

Ident. No.: 0020128734

~~€ 4,881,00~~ € 4.392,90

-10%

2+1
Year
Warranty*

ECO
energy-efficient
technology

Four integrated Watson-Marlow pumps



Discover our Bioprocessing Solutions HABITAT photo cell

- › Compact, space-saving design
- › Dimensions: 223 × 402 × 450 mm (W × D × H)
- › Latest technology sensors for better insights and more control
- › Large screen tablet for easy operation (10.4 inch)
- › Intuitive, easy-to-use control software that provides simple initial setup and intuitive user assistance throughout the process

Ident. No.: 0010007553

on request

2+1
Year
Warranty*

NEW!



Cooling and Shaking Incubator INC 125 FCS digital

- › High shaking speed from 10 to 300 rpm
- › Shaking stroke: 25 mm (optional: 20 mm) for optimum oxygen input
- › Decontamination mode (100 °C) for prevention of cross-contamination
- › Wide temperature range: 5 °C to 80 °C
- › Defrost function with continuous water drainage through the drain opening

Ident. No.: 0020135425 (SP20) | 0020136638 (SP25)

~~€ 7,798.00~~

€ 7.018,20

-10%

2+1
Year
Warranty*

ECO
energy-efficient
technology





Welcome to the IKA Knowledge Base

/// Find Insights and Hands-On Knowledge

Explore a collection of informative articles and practical examples showcasing a variety of applications in the laboratory field. As a result of our extensive collaboration with partners and customers from various industries, we have gained valuable know-how that we are excited to share with you. Our goal is to provide you with insights and hands-on knowledge that will help you make your processes more efficient.

Filter by Industry and Application to find the articles most valuable to you.



202603_OC22026_IWS_EN_EUR_94000336