

IKA

designed for scientists

EN



The Heart of Laboratory Innovation
Designed for You. Ready for the Future.

IKA Specials Q3 / 2026
Innovative Laboratory Solutions

imLab



www.imlab.eu - info@imlab.eu



+33(0)3 20 55 19 11



+32(0)16 73 55 72

Take advantage of our time-limited offers on essential lab equipment.

Get your discounts now!

The discounts are valid until September 30th, 2026.



NEW!



Electronic Pipette IKA PEN 5 – 100 µl

- › Available in three different volume ranges
- › Powered by a rechargeable lithium battery for operating a full working day on a single charge
- › With smart modes like reverse pipetting, multi-dispense or dilution of two liquids
- › With easy and straightforward self-calibration feature
- › Intuitive menu structure
- › Ergonomic design that fits naturally in your hand, reducing fatigue during repetitive tasks

2+1
Year
Warranty*

Ident. No.: 0020025102

~~€ 399~~ **€ 339.15**

-15%



IKA PEN App

NEW!

Fully manufactured
by robots "Made in
Germany"



Magnetic Stirrer RCT basic

- › Heating capacity of 800 W
- › IP Protection class IP 54
- › Compact design for less space required
- › Excellent magnetic coupling at up to 1500 rpm
- › With extensive multi-level safety concept for user and application
- › Functions expandable with USB dongle (available as accessory)

2+8
Year
Warranty*


Ident. No.: 0025008259

~~€ 898~~ **€ 763.30**

-15%*



*Also available as
anodized version



NEW!



*Also available as
anodized version

Magnetic Stirrer RET basic

- › Maximum speed of 1700 rpm
- › Easy connection to wireless sensors and thermometers
- › Integrated intelligent networking and digitalization features

Ident. No.: 0025008265

~~€ 1,298~~ € 1,103.30

-15%*

2+8
Year
Warranty*



NEW!

Fully functional
standalone pH meter



Complete Solutions
available at
www.ika.com

Combination with Wireless Sensor ETS-D7

- › For advanced temperature measurement and readability
- › Dual Control Interface with automatically synchronized values between magnetic stirrer and sensor
- › Rechargeable battery for 40 working hours (display 100% on), 80 working hours (display 50% on)

Ident. No.: 0025008370

~~€ 498~~ € 423.30

-15%

2+1
Year
Warranty*

NEW!



Overhead Stirrer EUROSTAR 60 digital

- › 23% smaller size than previous models
- › Torque Trend Measurement
- › Quick Sense Keyless Chuck
- › Cleanroom ISO 7/8 Classification
- › IP 54

Ident. No.: 0020105156

~~€ 1,648~~ € 1,400.80

-15%

2+1
Year
Warranty*

NEW!



Overhead Stirrer EUROSTAR 100 control

Additional features compared to digital version:

- › Viscosity measurement (mPas)
- › High precision torque measurement +/- 1 Ncm
- › Vibration sensor
- › Integrated Light
- › Left/Right rotation direction (interval mode)
- › Speed: Soft start in 3 modes
- › Chaotic Mixing Mode

Ident. No.: 0020105159

~~€ 2,998~~ € 2,548.30

-15%

2+1
Year
Warranty*

NEW!



Cooling and Shaking Incubator INC 125 FCS digital (SP25)

- › High shaking speed from 10 to 300 rpm
- › Shaking stroke: 25 mm (optional: 20 mm) for optimum oxygen input
- › Decontamination mode (100 °C) for prevention of cross-contamination
- › Wide temperature range: 5 °C to 80 °C
- › Defrost function with continuous water drainage through the drain opening

Ident. No.: 0020136638

~~€ 7,799~~ € 7,019.10

-10%

2+1
Year
Warranty*

ECO
energy-efficient
technology



With natural
refrigerant R290



Temperature Control for Viscometers CBC ROTAVISC SBS me-vi Package

- › Ready-to-connect solution for precise viscosity measurement of temperature-dependent samples
- › Heating and cooling bath covers a temperature range from -25 °C to 125 °C
- › Offers space for two standard beakers (600 ml)

Ident. No.: 0010014031

~~€ 7,568~~ € 6,811.20

-10%

2+1
Year
Warranty*

ECO
energy-efficient
technology



Magnetic Stirrer I-MAG Industry Stirrer

- › Powerful yet compact magnetic stirrer for applications in laboratory and production environments
- › Robust drive unit with IP protection class 64
- › Features magnetic coupling with high-performance neodymium magnets
- › Adjustable start-up speed that prevents the magnetic stirring bar from breaking off
- › In combination with the I-MAG SP set-up frames for large vessels up to 300 liters of water can be mixed efficiently

2+1
Year
Warranty*

Ident. No.: 0020014226

~~€ 4,201~~ **€ 3,780.90**

-10%



SP 150



SP 300

Accessories

I-MAG SP 150/SP 300 Platform extension

- › Stainless steel set-up frames for extension of footprint and load capacity of the I-MAG Industry stirrer
- › Designed for barrels up to approx. 150/300 liters and a maximum load of 200/400 kg

2+1
Year
Warranty*

Ident. No.: 0020102994 | 0020102996

~~€ 504 | € 756~~
€ 453.60 | € 680.40

-10%

Mix, heat, cool –
all in one



+ Free
PETTE 100 – 1000 µl



Thermoshaker MATRIX Orbital Delta F0.5

- › Large and clear display with a convenient menu structure
- › Aluminum die cast housing ensures secure positioning, excellent stability and everyday robustness
- › Quick and effective mixing of liquids with up to 2000 rpm and 3 mm shaker diameter
- › Programmable ramps allow for automatic execution of mixing and heating tasks

2+1
Year
Warranty*

Ident. No.: 0010006255

~~€ 2,736~~ **€ 2,462.40**

-10%