

An ergonomic feel combined with uncompromising precise performance: the IKA single fix and vario pipettes can be used for many applications in the modern laboratory. Colour coding simplifies the rapid selection of a suitable pipette. Resistance to impact, UV and chemicals ensures that the tip cone provides long-term reliable dosing. The IKA PETTE fix and vario are fully autoclavable, therefore ensuring easy sterilization. This guarantees fault-free, precise and safe operation.

IKA PETTE fix

/// Simple liquid handling for fixed volumes

ERGONOMIC DESIGN

IKA pipettes are delivered including exchangeable grips in different shapes and materials. The minimal frictional and spring forces also ensure dosing with ergonomic handling; without any effort.

DURABLE CONSTRUCTION

The tip cone and piston are coated with a high quality diamond-like carbon (DLC).

SIMPLE MAINTENANCE AND ADJUSTMENT

The volume unit with piston and seal can be disassembled and cleaned in just a few simple steps. If adjustment is necessary, it can be done without any additional tool.



- ① IKA PETTE fix 5 µl
EUR 144,00
Ident. No. 0020011219
- ② IKA PETTE fix 10 µl
EUR 144,00
Ident. No. 0020011220
- ③ IKA PETTE fix 20 µl
EUR 144,00
Ident. No. 0020011221
- ④ IKA PETTE fix 25 µl
EUR 144,00
Ident. No. 0020011222
- ⑤ IKA PETTE fix 50 µl
EUR 144,00
Ident. No. 0020011223
- ⑥ IKA PETTE fix 100 µl
EUR 144,00
Ident. No. 0020011224
- ⑦ IKA PETTE fix 200 µl
EUR 144,00
Ident. No. 0020011225
- ⑧ IKA PETTE fix 250 µl
EUR 144,00
Ident. No. 0020011226
- ⑨ IKA PETTE fix 500 µl
EUR 144,00
Ident. No. 0020011227
- ⑩ IKA PETTE fix 1000 µl
EUR 144,00
Ident. No. 0020011228



	IKA PETTE FIX	MAX. VOLUME TIP		IKA PETTE FIX	MAX. VOLUME TIP
①	5 µl	10 µl	⑥	100 µl	200 µl
②	10 µl	10 µl	⑦	200 µl	200 µl
③	20 µl	200 µl	⑧	250 µl	1000 µl
④	25 µl	200 µl	⑨	500 µl	1000 µl
⑤	50 µl	200 µl	⑩	1000 µl	1000 µl

① IKA PETTE vario 0,1 – 2 µl

EUR 219,00
Ident. No. 0020011210

② IKA PETTE vario 0,5 – 10 µl

EUR 219,00
Ident. No. 0020011211

③ IKA PETTE vario 2 – 20 µl

EUR 219,00
Ident. No. 0020011213

④ IKA PETTE vario 10 – 100 µl

EUR 219,00
Ident. No. 0020011214

⑤ IKA PETTE vario 20 – 200 µl

EUR 219,00
Ident. No. 0020011215

⑥ IKA PETTE vario 100 – 1000 µl

EUR 219,00
Ident. No. 0020011216

⑦ IKA PETTE vario 0,5 – 5 ml

EUR 219,00
Ident. No. 0020011217

⑧ IKA PETTE vario 1 – 10 ml

EUR 219,00
Ident. No. 0020011218



IKA PETTE VARIO MAX. VOLUME TIP

①	0,1 – 2 µl	10 µl
②	0,5 – 10 µl	10 µl
③	2 – 20 µl	200 µl
④	10 – 100 µl	200 µl

IKA PETTE VARIO MAX. VOLUME TIP

⑤	20 – 200 µl	200 µl
⑥	100 – 1000 µl	1000 µl
⑦	0,5 – 5 ml	5 ml
⑧	1 – 10 ml	10 ml



IKA LINEA



Pipette tips

Product	Description	Ident. No.	Price
ACCESSORIES			
LINEA	Pipette rack	0020018987	EUR 69,00
Tip xs bag	Pipette tip, 10 µl, transparent, 10 000 pcs.	0020017816	EUR 490,00
Tip xs tray	Pipette tip, 10 µl, transparent, 2304 pcs.	0020017818	EUR 156,00
Tip xs box +	Filter tip, 10 µl, sterile, transparent, 1920 pcs.	0020017820	EUR 300,00
Tip s bag	Pipette tip, 200 µl, transparent, 10 000 pcs.	0020017822	EUR 300,00
Tip s tray	Pipette tip, 200 µl, transparent, 1920 pcs.	0020017823	EUR 132,00
Tip s box +	Filter tip, 200 µl, sterile, transparent, 1920 pcs.	0020017824	EUR 300,00
Tip m bag	Pipette tip, 1000 µl, transparent, 5000 pcs.	0020017826	EUR 160,00
Tip m box	Pipette tip, 1000 µl, transparent, 1000 pcs.	0020017827	EUR 60,00
Tip m box +	Filter tip, 1000 µl, sterile, transparent, 1000 pcs.	0020017828	EUR 150,00
Tip l bag	Pipette tip, 5 ml, transparent, 2500 pcs.	0020017829	EUR 300,00
Tip l box	Pipette tip, 5 ml, transparent, 500 pcs.	0020017830	EUR 190,00
Tip xl bag	Pipette tip, 10 ml, transparent, 1000 pcs.	0020017831	EUR 325,00
Tip xl box	Pipette tip, 10 ml, transparent, 250 pcs.	0020017832	EUR 190,00
Pipette tip box xs/s	Stacking box, 6 pcs.	0020017834	EUR 90,00
Refill tray xs	Refill tray, 25 pcs.	0020017836	EUR 125,00
Refill tray s	Refill tray, 25 pcs.	0020017837	EUR 125,00
Pipette tip box l	Stacking box, 10 pcs.	0020017892	EUR 150,00



SOFTWARE

Pipette Calibration Software

Ident. No. 0020022141

IKA provides a user-friendly calibration software for planning, performing and recording all calibration tests.