


Traveler



Value and practicality for basic weighing

- Weigh-Thru draftshield for protection and stability
- Easily stackable storage
- Easy to install USB or RS232 interface options



- Application** Weighing (2 units)
- Display** LCD plus optional auxiliary display
- Operation** Mains (included), optional battery with programmable Auto Shut-Off
- Communication** Optional RS232 or USB interface with cable
- Construction** ABS plastic with stainless steel pan
- Design Features** Unique Weigh-thru draftshield designed with an easy to remove centre piece that provides protection from rough environments while maintaining weighing speed, designed to stack, shipping lock, adjustment lock, fixed integral weigh below hook



Model	Capacity	Readability (d)	Platform	Test Weight	Item N°	Price €
TA152	150 g	0.01 g	Ø 120 mm	150 g, M1, Included with balance	72202482	235
TA302	300 g	0.01 g	Ø 120 mm	300 g, M1, Included with balance	72202483	280
TA502	500 g	0.01 g	Ø 120 mm	300 g, M1	83997626	330
TA301	300 g	0.1 g	Ø 120 mm	300 g, M1, Included with balance	72202484	129
TA501	500 g	0.1 g	Ø 120 mm	300 g, M1	72202485	149
TA1501	1500 g	0.1 g	120 mm x 135 mm	1000 g, M1	72202486	220
TA3001	3000 g	0.1 g	120 mm x 135 mm	2000 g, M1	72202487	280
TA5000	5000 g	1 g	120 mm x 135 mm	3000 g, M1	72202489	149



Accessories for Navigator® and Traveler

Accessories	NV	NVT	TA	Item N°	Price €
RS232 interface kit	•	•	•	83032107	32
USB interface kit	•	•	•	83032108	32
Ethernet interface kit	•	•	•	83032109	67
Auxiliary Display TAD7			•	80251395	53
Battery rechargeable kit		•		83032106	40
Anti theft device	•	•	•	80850000	50
Carrying Case	•	•		30467963	on req.
Stacking Kit, x1	•	•		30467917	on req.
Stacking Kit, x4	•	•		30467918	on req.



FOR PRINTERS, CABLES AND SOFTWARE SEE P. 71-74

CS Series



Dependable results in a portable balance

- AC adapter included
- Large stainless steel pan included
- Auto Shut-off

Model	Capacity	Readability (d)	Platform	Test Weight	Item N°	Price €
CS200, 1 order = one pack of 6 pcs	200 g	0.1 g	145 mm x 133 mm	200 g, M1	83997617	360
CS600, 1 order = 1 pack of 6 pcs	600 g	0.1 g	145mm x 133mm	300g, M1	30422851	420
CS2000, 1 order = one pack of 6 pcs	2000 g	1 g	145 mm x 133 mm	1 kg, M1	83997618	360
CS5000, 1 order = one pack of 6 pcs	5000 g	1 g	145 mm x 133 mm	2 kg, M1	83997619	420



CSE Series



The perfect portable balance when cost matters and dependable results are a must!

- Large Weighing Platform
- Compact design
- Superior RF protection

Model	Capacity	Readability (d)	Platform	Test Weight	Item N°	Price €
CS200E, 1 order = 1 pack of 6 pcs	200 g	0.1 g	145 mm x 133 mm	200 g, M1	80500630	270
CS600E, 1 order = 1 pack of 6 pcs	600 g	0.1 g	145mm x 133mm	300g, M1	30422852	330
CS2000E, 1 order = 1 pack of 6 pcs	2000 g	1 g	145 mm x 133 mm	1 kg, M1	80500631	270
CS5000E, 1 order = 1 pack of 6 pcs	5000 g	1 g	145 mm x 133 mm	2 kg, M1	80500632	330



CL Series



Reliable portable scales when prices matters

- Storage made easy with a stackable design
- Integrated load cell lock
- Battery operation with auto shut-off (optional main adaptor)

Model	Capacity	Readability (d)	Platform	Item N°	Price €
CL201, 1 order = one pack of 6 pcs	200 g	0.1 g	Ø 120 mm	80010610	270
CL501, 1 order = one pack of 6 pcs	500 g	0.1 g	Ø 120 mm	80500640	330
CL2000, 1 order = one pack of 6 pcs	2000 g	1 g	Ø 120 mm	80010611	270
CL5000, 1 order = one pack of 6 pcs	5000 g	1 g	Ø 120 mm	80010612	330



Accessories for CS, CSE and CL Series

Accessories	CS	CSE	CL	Item N°	Price €
Hard shell carrying case	•	•	•	80010624	24
1050 ml bowl container	•	•	•	80850075	on req.
Stainless steel pan cover	•			80850076	on req.
Stainless steel pan cover. Reduces capacity by 80g			•	80010623	12
Power Adapter (Universal Plug Set Required)		•	•	46001724	25
Universal Plug Set		•	•	46001779	15

Precision Portable Balances



YA Gold Series



Reliable pocket-sized portable scale

- Unique hinged cover design
- High contrast backlit LCD
- Energy saving auto shut-off



Model	Capacity	Readability (d)	Platform	Item N°	Price €
YA102	100 g	0.01 g	50 mm x 70 mm	80251910	31
YA302	300 g	0.05 g	50 mm x 70 mm	80251911	27
YA501	500 g	0.1 g	50 mm x 70 mm	80251912	27

Emerald Series



Convenient, compact weighing capability in the palm of your hand

- Blue backlit luminescent display
- Protective snap-on cover



Model	Capacity	Readability (d)	Platform	Test Weight	Item N°	Price €
JE120, 1 order = one pack of 6 pcs	120 g	0.1 g	83 mm x 76 mm	100 g, M1	71148072	426
JE250, 1 order = one pack of 6 pcs	250 g	0.1 g	83 mm x 76 mm	200 g, M1	71148073	486
JE500, 1 order = one pack of 6 pcs	500 g	0.1 g	83 mm x 76 mm	200 g, + 100 g, M1	71148074	606

