

Ultra Turrax Homogenisator T 25 IKA - Promotie Package -



Best-seller IKA :
T25 digital



Silentstream



Dispergeer hulpstuk S 25 N 18 G



NEW // silentstream

without silentstream

- ⊖ formation of a whirl
- ⊖ introduction of air while dispersing
- ⊖ possible risk of running dry



with silentstream

- ✔ prevents the formation of a whirl
- ✔ less introduction of air into the medium
- ✔ prevents dry running risk



further advantages:

- ✔ dishwasher proofed
- ✔ sterilizable with all methods
- ✔ PEEK Material: resistance to high temperatures and chemicals

Toepassingen

- Homogenisatie van afvalwaterstalen
- Aanmaken van emulsies en dispersies
- Dispergeertaken onder vacuüm of druk
- Staalvoorbereiding voor medische diagnose

Eigenschappen

- Dispergeervolume : 1 tot 2.000 ml
- Digitaal toerentaldisplay
- Toerental regelbaar van 3.000 tot 25.000 tr/min
- Max. Viscositeitsbereik: 5.000 mPas



Designed to work perfectly

T 25 digital (ref. 0003725000)	1.955,00 €
+ Silentstream (ref. 0003754000)	+ 18,00 €
+ S25N-18G Dispergeer hulpstuk (ref. 0000593400)	+ 999,00 €
Catalogus Prijs :	= 2.954,00 €
Promotie Package (ref. 0009022900):	2.199,00 €

Algemene Verkoopsvoorwaarden:

Netto prijs in € / excl. BTW

Garantie: 2 jaar

Leveringstermijn: 1 à 3 dagen (voor Benelux & Frankrijk)

Aanbieding geldig zolang voorraad strekt