

# IKA

designed for scientists

## RW 20 digital Package

/// Data Sheet

The bestseller in the laboratory:

- With digital display
- Robust, slimline, ergonomic design
- Technical improvements on the trusted RW 20 series designs
- With constant power-drive
- Two speed ranges for universal use from 60 - 2.000 rpm
- Push-through agitator shafts (only when stationary)
- For stirring quantities of up to 20 l (H<sub>2</sub>O)



RW 20 digital

Deux plages de vitesse pour une utilisation universelle de 60 à 2.000 min<sup>-1</sup>

+



R 1826

Statif sur pied

+



R 1342

Tige à hélice

+



R 182

Noix de serrage

RW 20 digital Package

Prix catalogue ~~€ 1.689,00~~

Prix de l'ensemble € 1.299,00

[www.imlab.eu](http://www.imlab.eu)

[info@imlab.be](mailto:info@imlab.be) - Tél. : +32 (0)16 73 55 72

**imLab**

[www.imlab.eu](http://www.imlab.eu)

[info@imlab.fr](mailto:info@imlab.fr) - Tél. : +33 (0)3 20 55 19 11



designed for scientists

## Technical Data

Stirring quantity max. per stirring position (H2O) [l]	20
Motor rating input [W]	70
Motor rating output [W]	35
Speed display	LED
Speed range [rpm]	60 - 2000
Viscosity max. [mPas]	10000
Output max. at stirring shaft [W]	26
Permissible ON time [%]	100
Torque max. at stirring shaft [Ncm]	150
Torque max. at stirring shaft at 60 1/min (overload) [Ncm]	300
Torque max. at stirring shaft at 100 1/min [Ncm]	150
Torque max. at stirring shaft at 1.000 1/min [Ncm]	24
Speed range I (50 Hz) [rpm]	60 - 500
Speed range II (50 Hz) [rpm]	240 - 2000
Speed range I (60 Hz) [rpm]	72 - 600
Speed range II (60 Hz) [rpm]	288 - 2400
Speed control	stepless
Setting accuracy speed [ $\pm$ rpm]	1
deviation of speed measurement [ $\pm$ rpm]	30
Stirring element fastening	chuck
Chuck range diameter [mm]	0.5 - 10
Hollow shaft, inner diameter [mm]	10.5
Hollow shaft (push-through - when stopped)	yes
Fastening on stand	extension arm
Extension arm diameter [mm]	13
Extension arm length [mm]	160
Speed control	mechanical
Nominal torque [Nm]	1.5
Dimensions (W x H x D) [mm]	88 x 294 x 212
Weight [kg]	3.1
Permissible ambient temperature [°C]	5 - 40
Permissible relative humidity [%]	80
Protection class according to DIN EN 60529	IP 20
Voltage [V]	220 - 240 / 100 - 115
Frequency [Hz]	50/60
Power input [W]	72

[www.imlab.eu](http://www.imlab.eu)

info@imlab.be - Tél. : +32 (0)16 73 55 72

**imLab**

[www.imlab.eu](http://www.imlab.eu)

info@imlab.fr - Tél. : +33 (0)3 20 55 19 11

Imlab sprl - Oude Vijvers, 1  
3370 Boutersem - Belgique

Imlab sarl - Centre d'Affaires de l'Horlogerie  
48 rue des Canoniers - 59000 Lille - France