

## C-MAG HS 7 Package

**List price:**  
1.187,00 €

**Package price:**  
**849,00 €**

Ident. No. 9015900

C-MAG HS 7  
Magnetic stirrer with heating  
new: with ceramic plate



- **Fixed safety circuit**  
Shuts off heater if temperature rises above 550 °C
- **Hot top indicator**  
Light flashes to warn user that the surface is hot to the touch
- **Digital display**  
Allows for accurate and repeatable temperature settings
- **Powerful motor**  
Stirs viscous materials with a maximum volume of 15 l (H<sub>2</sub>O)
- **Strong magnet**  
High magnetic adhesion prevents stir bar decoupling
- **Microprocessor control**  
Sends feedback to the motor and hotplate to maintain constant speed and set temperature
- **Ceramic top**  
Seamless one-piece glass-ceramic top provides excellent resistance to chemicals
- **Temperature controller**  
Connection for contact thermometer enables precise temperature control of the medium
- **Elevated control panel**  
Minimizes contact with spills and allows for easy reading and adjustment of controls



- + ETS-D5 Electronic contact thermometer
- + H 16 V Support rod
- + H 38 Holding rod
- + H 44 CLAMP